

Selecting and ordering earmolds



signia

Tips for venting selection	2
Receiver strength recommendations based on general audiogram	2
Material selection.....	3
Acrylic.....	3
Silicone	3
Style selection	3
Canal type selection for acrylic RIC 3.0 earmolds.....	4
Canal lengths defined	5
Good impressions lead to good fittings:	5
BTE tubing options.....	5
Standard tubing in single bend or double bend.....	5
Dry tube	6
Quick-connect tubing system	6
Tube lock.....	6
ThinTube 3.0.....	6



signia

Tips for venting selection



To maximize hearing aid processing benefit, please use the most comfortable closed coupling method available

General guidelines that can be followed if closed coupling fitting is not tolerated by the client

- Normal hearing below 1kHz: 3.0mm vent is recommended
- Mild hearing loss below 1kHz: 2.0-2.5mm vent is recommended
- Moderate hearing loss below 1kHz: 1.3-1.6mm vent is recommended
- Severe hearing loss below 1kHz: pressure or no vent is recommended
- Profound hearing loss below 1kHz: no vent/ closed fitting is recommended

Receiver strength recommendations based on general audiogram

- S receiver (110/46) covers hearing loss up to 50dB in the low frequencies and 80dB in the high frequencies*
 - M receiver (119/60) covers hearing loss up to 70dB in the low frequencies and 95dB in the high frequencies*
 - P receiver (122/65) covers hearing loss up to 90dB in the low frequencies and 105dB in the high frequencies*
 - HP receiver (131/75) covers hearing loss up to 100dB in the low frequencies and 115dB in the high frequencies*
- The HP form factor is a unique receiver type that is built with a custom shell and faceplate in flesh-tone colors.



*Assumes sensorineural hearing loss. More gain may be needed when a conductive component is present



signia

Material selection

Acrylic

- Ideal for seniors as the outer ear becomes softer with age. Harder materials are less likely to pull on softer skin
- Durable, long lasting and not apt to discolor and change over time
- Easy to clean
- Easy to modify in the office

Silicone

- Ideal for young children with growing ears and active patients
- May offer better retention

Style selection



RIC

The RIC and RIC Lock are the most cosmetically appealing options.



Canal Foil

The canal models, canal, canal foil, canal foil lock are built larger in the concha area than the RIC models



Skeleton

Skeleton earmolds have a hollow center but sit in the concha bowl like full shell. These molds have a more cosmetic-friendly appearance and often are more comfortable.



RIC Lock

The RIC and RIC Lock are the most cosmetically appealing options. The difference is that the RIC lock comes with a canal lock.



Canal Foil Lock

Anytime the word "Foil" is used, that means the inside of the mold is hollowed out and the earmold will accommodate receiver strengths for which it was not built. This is great for users on the borderline of fitting ranges, or for those whose hearing loss will progress. The canal models, canal, canal foil, canal foil lock, canal lock are built larger in the concha area than the RIC models.



Skeleton Foil

Anytime the word "Foil" is used, that means the inside of the mold is hollowed out and the earmold will accommodate receiver strengths for which it was not built. This is great for users on the borderline of fitting ranges, or for those whose hearing loss will progress.



Foil RIC

Anytime the word "Foil" is used, that means the inside of the mold is hollowed out and the earmold will accommodate receiver strengths for which it was not built. This is great for users on the borderline of fitting ranges, or for those whose hearing loss will progress.



Half Shell

Half-shell earmolds cover only half of the concha bowl.



3/4 Skeleton

For 3/4 Skeleton the rim of material in the concha bowl area is reduced to 3/4 lock.



RIC Foil Lock

Anytime the word "Foil" is used, that means the inside of the mold is hollowed out and the earmold will accommodate receiver strengths for which it was not built. This is great for users on the borderline of fitting ranges, or for those whose hearing loss will progress. The canal lock is built larger in the concha area than the RIC models.



Helix Lock

Helix locks have a small extension out of the canal along the top of the concha to the helix for better retention.



1/2 Skeleton

For 1/2 Skeleton the rim of material in the concha bowl area is reduced even further to 1/2 lock.



Canal

The canal models, canal, canal foil, canal foil lock are built larger in the concha area than the RIC models.



Half Shell Foil

Anytime the word "Foil" is used, that means the inside of the mold is hollowed out and the earmold will accommodate receiver strengths for which it was not built. This is great for users on the borderline of fitting ranges, or for those whose hearing loss will progress.



Semi Skeleton

The semi skeleton earmold removes part of the outer lock and leaves the final mold with a lower canal and helix lock.



Canal Lock

The canal models, canal, canal foil, canal foil lock, canal lock are built larger in the concha area than the RIC models.



Full Shell

Full-shell earmolds generally cover the entire concha bowl and provide a nice acoustic seal.

Note: Earmolds will ship with RIC when ordered together.



signia

Canal type selection for acrylic RIC 3.0 earmolds

Acrylic RIC 3.0 earmolds have two choices for canal type, long or short. These choices determine the canal lengths available as well as the wax guard configuration.

Long

1 | Select Long Canal Type

Canal Type *

Long

Please select an option

Long

Short

2 | Select Wax Guard preference for Earmold (3 options)

Choosing the NanoCare Wax Guard ensures that both the miniReceiver 3.0 and RIC Earmold use the same wax guards. Ease of use for you and your patients.

Waxguard *

NanoCare Wax Guard

Please select an option

NanoCare Wax Guard

No Wax Guard

Quickguard Wax Guard

*Auto default is NanoCare Wax Guard.
The miniReceiver 3.0 uses the NanoCare Wax Guard.*

3 | Select Canal Length (5 options)

Canal length *

Medium

Please select an option

Customer specified

Deep

Long

Medium

Short

Short

1 | Select Short Canal Type

Canal Type *

Long

Please select an option

Long

Short

2 | Select Wax Guard Option (1 option)

Waxguard *

NanoCare Wax Guard

Please select an option

NanoCare Wax Guard

*NanoCare Wax Guard is the only option for Short Canal type.
The miniReceiver 3.0 uses the NanoCare Wax Guard.*

3 | Select Canal Length (3 options)

Canal length *

Medium

Please select an option

Customer specified


Medium

Short



signia

Canal lengths defined



Good impressions lead to good fittings:

Complete impression, at least 2mm past the second bend and free of voids

Consider ear anatomy and size for earmold and receiver strength options

BTE tubing options

Standard tubing in single bend or double bend

Single Bend



Double Bend



Standard tubing size is indicated by the outside diameter of the tubing:

3.1 = standard size. For mild or moderate hearing loss

3.3 = thick tubing. For severe to profound hearing loss

3.6 = heavy walled tubing. For profound hearing loss

Dry tube

A dry tube option is available in all diameters, thicknesses and single or double bend options. The dry tube is useful for clients who present with moisture accumulation in the standard tubing. Dry tube is porous, yet acoustically transparent material, that allows moisture to evaporate naturally.


Quick-connect tubing system

The quick-connect tubing system requires no glue, so it is not necessary to ream the earmold when replacing the tubing. This system also allows the tubing to turn in the earmold providing a more natural tubing position for your client. The quick-connect tubing system is offered for all standard BTE earmolds, both acrylic and silicone.



Tube lock

A tube lock tube is default on all silicone earmolds ordered with standard tubing. This system eliminates the need for glue, making the replacement of tubing quick and easy. The lock mechanism is a clear plastic ring which does not discolor with age.



ThinTube 3.0

The improved bending of the tubing and outlet preserves the acoustic coupling in the ear canal, prevents lateral migration, and reduces the need for a concha lock.