




Signia Charger-Übersicht



	Multi Charger 21026076	Silk Mobile Charger 21025470	Insio Charger 21025145	Active IX Mobile Charger 21028352	Styletto IX Mobile Charger 21031113	Pure Mobile Charger 10993119	Pure Charger 10981896	Motion P Charger 10981894	Motion SP Charger 10981895	Styletto IX Dry&Clean Charger 21007594	Pure Dry&Clean Charger 10991651	Motion P Dry&Clean Charger 10991639	Motion SP Dry&Clean Charger 10991645	Induktiv-Charger II* 10966943
Silk Charge&Go	-	●	-	-	-	-	-	-	-	-	-	-	-	-
Pure Charge&Go (T)	●	-	-	-	-	●	●	-	-	-	●	-	-	-
Styletto IX	●	-	-	-	●	-	-	-	-	●	-	-	-	-
Insio Charge&Go	-	-	●	-	-	-	-	-	-	-	-	-	-	-
Active Pro IX/ Active IX	-	-	-	●	-	-	-	-	-	-	-	-	-	-
Motion Charge&Go	-	-	-	-	-	-	-	-	-	-	-	-	-	●
Motion Charge&Go P	●	-	-	-	-	-	-	●	-	-	-	●	-	-
Motion Charge&Go SP	●	-	-	-	-	-	-	-	●	-	-	-	●	-

● voll kompatibel

*Alternativ: Induktiv-Charger E (schwarz/offen) mit der Sachnummer 10991800