

 EarWear 3.0 – Discreto. Seguro. Cómodo.

Tabla de compatibilidad de EarWear

	ThinTube 3.0	miniReceiver 3.0	Eartips 3.0 y Sleeves 3.0	Earmold 3.0	miniReceiver 2.0	Click Domes y Click Sleeves	Click Mold
Serie Motion IX	● ***	–	●	●	–	–	–
Serie Pure IX	–	●	●	●	–	–	–
Serie Styletto IX	–	●	●	●	–	–	–
Serie Silk IX	–	–	● *	● **	–	–	–
Serie Active Pro IX	–	–	● *	–	–	–	–
Serie Silk X	–	–	–	–	–	●	–

* Compatible solo con Sleeve 3.0

** Los earmolds personalizados para la serie Silk IX no son intercambiables con el Earmold 3.0 estándar

*** Motion Charge&Go SP IX y Motion Charge&Go SP X no son compatibles con ThinTube 3.0