

EarWear compatibility table

	ThinTube 3.0	miniReceiver 3.0	Eartips 3.0 & Sleeves 3.0	Earmold 3.0	miniReceiver 2.0	Click Domes & Click Sleeves	Click Mold
Motion IX series	● ***	–	●	●	–	–	–
Pure IX series	–	●	●	●	–	–	–
Styletto IX series	–	●	●	●	–	–	–
Silk IX series	–	–	● *	● **	–	–	–
Active Pro IX series	–	–	● *	–	–	–	–
Styletto AX series	–	●	●	●	–	–	–
Pure AX series	–	●	●	●	–	–	–
Motion X series	● ***	–	●	●	–	–	–
Pure 312 X	–	–	–	–	●	●	●
Silk X series	–	–	–	–	–	●	–

* Compatible with Sleeve 3.0 only
** Custom earmolds for Silk IX series are not interchangeable with standard Earmold 3.0
*** Motion Charge&Go SP IX and Motion Charge&Go SP X are not compatible with ThinTube 3.0