



Signia Charger Compatibility

Hearing aid model	Charger model			
	Portable Chargers	Desktop Chargers	Dry&Clean Chargers	Multi Charger Compatibility
Integrated Xperience				
Pure Charge&Go BCT IX (and CROS)	Pure Portable Charger	Desktop Charger Pure	Dry&Clean Charger Pure	Yes
Pure Charge&Go T IX	Pure Portable Charger	Desktop Charger Pure	Dry&Clean Charger Pure	Yes
Pure Charge&Go IX (and CROS)	Pure Portable Charger	Desktop Charger Pure	Dry&Clean Charger Pure	Yes
Styletto IX (and CROS)	Styletto IX Charger	-	Dry&Clean Charger Styletto	-
Motion Charge&Go M IX	-	Desktop Charger Motion M	-	Yes
Motion Charge&Go P IX	-	Desktop Charger Motion P	Dry&Clean Charger Motion P	Yes
Motion Charge&Go SP IX	-	Desktop Charger Motion SP	Dry&Clean Charger Motion SP	Yes
Active Pro IX (and CROS)	Active IX Charger	-	-	-
Silk Charge&Go IX (and CROS)	Silk Charger	-	-	-
Insio Charge&Go CIC IX	Portable Charger Insio CIC	-	-	-
Augmented Xperience				
Pure Charge&Go AX (and CROS)	Pure Portable Charger	Desktop Charger Pure	Dry&Clean Charger Pure	Yes
Pure Charge&Go T AX	Pure Portable Charger	Desktop Charger Pure	Dry&Clean Charger Pure	Yes
Styletto AX (and CROS)	Styletto AX Charger	Desktop Charger Pure	Dry&Clean Charger Styletto	-
Insio Charge&Go AX	-	Desktop Charger Insio AX	-	-
Xperience				
Motion Charge&Go X	-	Inductive Charger II, Inductive Charger E* *for 2X/1X	-	-
Motion Charge&Go P X	-	Desktop Charger Motion P	Dry&Clean Charger Motion P	Yes
Motion Charge&Go SP X	-	Desktop Charger Motion SP	Dry&Clean Charger Motion SP	Yes