

Overzicht per hoortoestel Tubes & Receivers Tips & Sleeves



Pure C&G BCT IX

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Pure C&G AX

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Active IX & Pro IX

- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Pure C&G T AX

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Styletto IX

- MiniReceiver 3.0 (S,M,P)
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Pure 312 AX

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Pure C&G IX

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Styletto AX

- MiniReceiver 3.0 (S,M,P)
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Pure C&G T IX

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Motion C&G P X

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Silk C&G IX

- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)

Overzicht per hoortoestel Tubes & Receivers Tips & Sleeves



Motion C&G X

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Motion C&G NX

- ThinTube 2.0 met koppelstukje
- Tips (Open, Tulip, Closed, Double)



Silk X

- ClickSleeves (Vented, Closed)



Pure C&G NX

- MiniReceiver 2.0 (S,M,P)
- MiniReceiver 2.0 HP
- ClickDomes (Open, Tulip, Closed, Double)
- ClickSleeves (Vented, Closed)



Styletto X

- MiniReceiver 2.0 (S,M,P)
- ClickDomes (Open, Tulip, Closed, Double)
- ClickSleeves (Vented, Closed)



Pure 312 NX

- MiniReceiver 2.0 (S,M,P)
- MiniReceiver 2.0 HP
- ClickDomes (Open, Tulip, Closed, Double)
- ClickSleeves (Vented, Closed)



Pure C&G X

- MiniReceiver 2.0 (S,M,P)
- MiniReceiver 2.0 HP
- ClickDomes (Open, Tulip, Closed, Double)
- ClickSleeves (Vented, Closed)



Pure 10 NX

- MiniReceiver 2.0 (S,M,P)
- ClickDomes (Open, Tulip, Closed, Double)
- ClickSleeves (Vented, Closed)



Pure 312 X

- MiniReceiver 2.0 (S,M,P)
- MiniReceiver 2.0 HP
- ClickDomes (Open, Tulip, Closed, Double)
- ClickSleeves (Vented, Closed)



Motion NX

- ThinTube 2.0 met koppelstukje
- Tips (Open, Tulip, Closed, Double)

Overzicht kleurenmogelijkheden per toestel



Pure C&G BCT IX - Pure C&G IX - Pure C&G T IX - Pure C&G AX - Pure C&G T AX - Pure 312 AX



Active IX & Active Pro IX




Styletto IX




Styletto X



Styletto AX



Silk IX – Silk X



Pure C&G X - Pure 312 X - Motion C&G X - Motion C&G P X - Motion C&G SP X



Motion C&G NX



Pure C&G NX - Pure 312 NX - Pure 10 NX - Motion NX



Overzicht filters bij de juiste Tubes & Receivers IHO & Molds

Voor itemnummers klik hier of ga naar pagina 9



WaxGuard 3.0 NanoCare

- ThinTube 3.0 (Standaard)
- ThinTube 3.0 (Power)
- MiniReceiver 3.0 (S,M,P)
- Extended Earmold 3.0
- Silk C&G IX
- Active IX & Active Pro IX

Voor instructie scan of klik op de QR code



ThinTube 3.0

WaxGuard MiniReceiver

- MiniReceiver 2.0
- Silk X, Silk NX




MiniReceiver 2.0

QuickGuard

- ClickMold
- MiniReceiver 2.0 HP
- Insio NX
- Insio C&G AX
- Insio IX IIC/CIC

Voor instructie scan of klik op de QR code



Clickmold

WaxGuard 3.0 Cerustop XL

- MiniReceiver 3.0 HP



MiniReceiver 3.0 HP

MF1 Microfoonfilter

- Insio IX IIC/CIC
- Insio C&G AX
- Insio NX CIC model
- in Zwart, Donkerbruin, Bruin, Mokka en Beige.
(Insio NX CIC alleen in Bruin en Mokka)
- MF1 Removal Tool ook los te bestellen

Voor instructie scan of klik op de QR code



Insio C&G AX

Overzicht tools


bij de juiste Tubes & Receivers IHO & Molds

Voor itemnummers klik hier of ga naar pagina 9



Inσιο IX Quick-Check card

- Afdruk meten of deze geschikt is voor productie Inσιο IX
Voor instructie scan of klik op de QR code




Meettool Earwear 3.0

- Lengte Tube 3.0 en miniReceiver 3.0 meten




Connexx EarMold 3.0 Measuring card

- Afdruk meten of deze geschikt is voor Maatwerk
t.b.v. Earwear 3.0




Earmold 3.0 Removal Tool

- Verwijderen Earmold 3.0 (kort)
Voor instructie scan of klik op de QR code



Earmold 3.0 Inserting and Removal Tool

- Plaatsen en verwijderen van Tube in Earmold 3.0 (kort en lang)
en miniReceiver 3.0
Voor instructie scan of klik op de QR code




ClickMold Removal Tool

- Verwijderen Clickmold van een miniReceiver 2.0



MF1 Removal Tool

- Verwijderen en plaatsen MF1 Filter
Voor instructie scan of klik op de QR code



Tool Removal H2770 BLU

- Verwijderen van de toonbocht bij Motion X hoortoestellen
- Verwijderen van de miniReceiver 2.0



Overzicht Conchalock 3.0 bij de juiste Tubes

Voor itemnummers klik hier of ga naar pagina 9



ConchaLock 3.0 universeel


- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- MiniReceiver 3.0 (S,M,P)

Voor instructie scan of klik op de QR code



Conchalock 3.0 maat S

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- MiniReceiver 3.0 S



Conchalock 3.0 maat M

- MiniReceiver 3.0 M



Conchalock 3.0 maat P

- MiniReceiver 3.0 P



Overzicht


Tips & Sleeves bij de juiste Tubes en Receivers

Voor itemnummers klik hier of ga naar pagina 9



EarTips 3.0 (Open, Tulip)

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- MiniReceiver 3.0 (S,M,P)
- Silk C&G IX
- Active IX & Active Pro IX



Sleeves 3.0 XS,S,M,L (Vented, Closed, Power)

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- MiniReceiver 3.0 (S,M,P)
- Silk C&G IX
- Active IX & Active Pro IX



ClickDomes 2.0 (Open, Tulip, Closed, Double)

- MiniReceiver 2.0 (S,M,P)




ClickSleeves 2.0 XS,S,M,L (Vented, Closed)

- MiniReceiver 2.0 (S,M,P)
- Silk X, Silk NX



Tips (Open, Tulip, Closed, Double)

- Life Tube
- ThinTube 2.0 met koppelstukje



Overzicht

ThinTubes en miniReceivers

Voor itemnummers klik hier of ga naar pagina 10




BE

MiniReceivers 3.0

- S Lengte 0 t/m 5
- M Lengte 0 t/m 5
- P Lengte 0 t/m 5
- HP Lengte 0 t/m 5 (via maatwerk bestellen)

Voor instructie scan of klik op de QR code



S

M

P

HP

ThinTubes 3.0

- Standaard 0,9 mm - Lengte 0 t/m 4
- Power 1,4 mm - Lengte 0 t/m 4



ThinTube 3.0
standaard 0,9 mm



ThinTube 3.0
Power 1,4 mm

MiniReceivers 2.0

- S Lengte 0 t/m 3
- M Lengte 0 t/m 4
- P Lengte 0 t/m 3
- HP Lengte 0 t/m 3 (via maatwerk bestellen)



S

M

P

HP

ThinTubes

- ThinTube 2.0 met koppelstukje
Lengte 1 t/m 4
- Life Tube
Lengte 1 t/m 6



ThinTube met koppelstukje



LifeTube

Overzicht filters

WaxGuard 3.0 NanoCare	10993649	WAXGUARD 3.0 NANOCARE CXX
WaxGuard MiniReceiver	10962332	WAXGUARD MINIRECEIVER (CXX)
QuickGuard	10962504	QUICKGUARD PACK
WaxGuard 3.0 XL NanoCare	10999334	WAXGUARD 3.0 CERUSTOP XL CXX G
MF1 Microfoonfilter	10947542	SERVICE SET (MF1) MOC
	10947776	SERVICE SET (MF1) BRN
	21016200	SERVICE SET (MF1) BLK
	21016201	SERVICE SET (MF1) SK
	21016202	SERVICE SET (MF1) DBN

Overzicht tools

Meettool Earwear 3.0	10991403	GAUGE TEMPLATE (PR) MINIREC 3.0
Connexx EarMold 3.0 Measuring card	21007585	EARMOLD 3.0 MEASURING CARD
Earmold 3.0 Removal Tool	10969958	TOOL HANDLING 3.0
Earmold 3.0 Inserting and Removal Tool	10967340	TOOL REMOVAL 3.0
ClickMold Removal Tool	10598984	SET CLICK MOLD SERVICE
MF1 Removal Tool	10147702	TOOL WECHSELSTIFT MF1
Tool Removal H2770 BLU	10600047	TOOL REMOVAL H2770 BLU

Overzicht Conchalock 3.0

ConchaLock 3.0 universeel	10993677	PACK CONCHA LOCK 3.0
Conchalock 3.0 maat S	21024885	PACK CONCHA LOCK SLEEVE 3.0 S
Conchalock 3.0 maat M	21024886	PACK CONCHA LOCK SLEEVE 3.0 M
Conchalock 3.0 maat P	21024887	PACK CONCHA LOCK SLEEVE 3.0 P

Overzicht Tips & Sleeves

EarTips 3.0 (Open, Tulip)	10993678	EARTIP 3.0 5MM OPEN CXX
	10993680	EARTIP 3.0 7MM OPEN CXX
	10993681	EARTIP 3.0 10MM OPEN CXX
	10993682	EARTIP 3.0 8MM TULIP CXX
	10993683	EARTIP 3.0 12MM TULIP CXX
Sleeves 3.0 XS,S,M,L (Vented, Closed, Power)	10992863	SLEEVE 3.0 VENTED XS CXX
	10992864	SLEEVE 3.0 VENTED S CXX
	10992865	SLEEVE 3.0 VENTED M CXX
	10992866	SLEEVE 3.0 VENTED L CXX
	10992867	SLEEVE 3.0 CLOSED XS CXX
	10992868	SLEEVE 3.0 CLOSED S CXX
	10992869	SLEEVE 3.0 CLOSED M CXX
	10992870	SLEEVE 3.0 CLOSED L CXX
	10997192	SLEEVE 3.0 POWER XS CXX
	10997193	SLEEVE 3.0 POWER S CXX
	10997194	SLEEVE 3.0 POWER M CXX
10997195	SLEEVE 3.0 POWER L CXX	
ClickDomes 2.0 (Open, Tulip, Closed, Double)	10426012	CLICK DOME 8MM CLOSED CXX TR
	10426013	CLICK DOME 10MM CLOSED CXX TR
	10426014	CLICK DOME 8MM OPEN CXX TR
	10426015	CLICK DOME 10MM OPEN CXX TR
	10426016	CLICK DOME 4MM OPEN CXX TR
	10426017	CLICK DOME 6MM CLOSED CXX TR
	10426018	CLICK DOME 8MM DOUBLE CXX TR
	10426019	CLICK DOME 10MM DOUBLE CXX TR
	10670790	CLICK DOME SEMI OPEN CXX TR
ClickSleeves 2.0 XS,S,M,L (Vented, Closed)	10951146	CLICK SLEEVE 2.0 S VENTED CXX
	10951147	CLICK SLEEVE 2.0 S CLOSED CXX
	10951148	CLICK SLEEVE 2.0 M VENTED CXX
	10951149	CLICK SLEEVE 2.0 M CLOSED CXX
	10951150	CLICK SLEEVE 2.0 L VENTED CXX
	10951151	CLICK SLEEVE 2.0 L CLOSED CXX
	10951152	CLICK SLEEVE 2.0 XS VENTED CXX
	10951153	CLICK SLEEVE 2.0 XS CLOSED CXX
Tips (Open, Tulip, Closed, Double)	10169328	LIFE TIP REFILL 4 MM
	10169329	LIFE TIP REFILL 6MM
	10054916	LIFE TIP REFILL 8MM
	10169330	LIFE TIP REFILL 8MM FULL
	10054917	LIFE TIP REFILL 10MM
	10175566	LIFE TIP DOUBLE REFILL 8/10MM
	10175715	LIFE TIP DOUBLE REFILL 10/12MM
	10962754	LIFE TIP SEMI OPEN REFILL

Overzicht ThinTubes en miniReceivers

MiniReceivers 3.0	10992548	MINIRECEIVER 3.0 0R S	MiniReceivers 2.0	10669907	MINIRECEIVER 2.0 MT. 0LS	
	10992549	MINIRECEIVER 3.0 0L S		10669908	MINIRECEIVER 2.0 MT. 0RS	
	10992550	MINIRECEIVER 3.0 1R S		10669909	MINIRECEIVER 2.0 MT. 1LS	
	10992552	MINIRECEIVER 3.0 1L S		10669910	MINIRECEIVER 2.0 MT. 1RS	
	10992553	MINIRECEIVER 3.0 2R S		10669911	MINIRECEIVER 2.0 MT. 2LS	
	10992554	MINIRECEIVER 3.0 2L S		10669912	MINIRECEIVER 2.0 MT. 2RS	
	10992555	MINIRECEIVER 3.0 3R S		10669913	MINIRECEIVER 2.0 MT. 3LS	
	10992556	MINIRECEIVER 3.0 3L S		10669914	MINIRECEIVER 2.0 MT. 3RS	
	10992557	MINIRECEIVER 3.0 4R S		10824820	MINIRECEIVER 2.0 MT. 0LM	
	10992558	MINIRECEIVER 3.0 4L S		10824821	MINIRECEIVER 2.0 MT. 0RM	
	21023207	MINIRECEIVER 3.0 5R S		10824822	MINIRECEIVER 2.0 MT. 1LM	
	21023211	MINIRECEIVER 3.0 5L S		10824823	MINIRECEIVER 2.0 MT. 1RM	
	10992559	MINIRECEIVER 3.0 0R M		10824824	MINIRECEIVER 2.0 MT. 2LM	
	10992560	MINIRECEIVER 3.0 0L M		10824825	MINIRECEIVER 2.0 MT. 2RM	
	10992561	MINIRECEIVER 3.0 1R M		10824826	MINIRECEIVER 2.0 MT. 3LM	
	10992562	MINIRECEIVER 3.0 1L M		10824827	MINIRECEIVER 2.0 MT. 3RM	
	10992563	MINIRECEIVER 3.0 2R M		10933054	MINIRECEIVER 2.0 MT. 4LM	
	10992564	MINIRECEIVER 3.0 2L M		10933053	MINIRECEIVER 2.0 MT. 4RM	
	10992565	MINIRECEIVER 3.0 3R M		10669923	MINIRECEIVER 2.0 MT. 0LP	
	10992566	MINIRECEIVER 3.0 3L M		10669924	MINIRECEIVER 2.0 MT. 0RP	
	10992567	MINIRECEIVER 3.0 4R M		10669925	MINIRECEIVER 2.0 MT. 1LP	
	10992568	MINIRECEIVER 3.0 4L M		10669926	MINIRECEIVER 2.0 MT. 1RP	
	21023215	MINIRECEIVER 3.0 5R M		10669927	MINIRECEIVER 2.0 MT. 2LP	
	21023216	MINIRECEIVER 3.0 5L M		10669928	MINIRECEIVER 2.0 MT. 2RP	
	10992569	MINIRECEIVER 3.0 0R P		10669929	MINIRECEIVER 2.0 MT. 3LP	
	10992570	MINIRECEIVER 3.0 0L P		10669930	MINIRECEIVER 2.0 MT. 3RP	
	10992571	MINIRECEIVER 3.0 1R P		ThinTubes	10828291	CXX TT SIZE 1L
	10992572	MINIRECEIVER 3.0 1L P			10828292	CXX TT SIZE 1R
	10992573	MINIRECEIVER 3.0 2R P			10828293	CXX TT SIZE 2L
	10992574	MINIRECEIVER 3.0 2L P			10828294	CXX TT SIZE 2R
	10992575	MINIRECEIVER 3.0 3R P			10828295	CXX TT SIZE 3L
	10992576	MINIRECEIVER 3.0 3L P			10828296	CXX TT SIZE 3R
	10992577	MINIRECEIVER 3.0 4R P			10828297	CXX TT SIZE 4L
	10992578	MINIRECEIVER 3.0 4L P			10828298	CXX TT SIZE 4R
	21023217	MINIRECEIVER 3.0 5R P			10054892	LIFE TUBE REFILL 1 BLUE
	21023218	MINIRECEIVER 3.0 5L P			10054891	LIFE TUBE REFILL 1 RED
	ThinTubes 3.0	10991609			Set ThinTube 3.0 S0	10054894
10991610		Set ThinTube 3.0 S1	10054893		LIFE TUBE REFILL 2 RED	
10991611		Set ThinTube 3.0 S2	10054896		LIFE TUBE REFILL 3 BLUE	
10991612		Set ThinTube 3.0 S3	10054895		LIFE TUBE REFILL 3 RED	
10991613		Set ThinTube 3.0 S4	10054898		LIFE TUBE REFILL 4 BLUE	
10992492		Set ThinTube 3.0 S0 P	10054897		LIFE TUBE REFILL 4 RED	
10992493		Set ThinTube 3.0 S1 P	10054900		LIFE TUBE REFILL 5 BLUE	
10992494		Set ThinTube 3.0 S2 P	10054899		LIFE TUBE REFILL 5 RED	
10992495		Set ThinTube 3.0 S3 P	10054902	LIFE TUBE REFILL 6 BLUE		
10992496		Set ThinTube 3.0 S4 P	10054901	LIFE TUBE REFILL 6 RED		