

# Signia CROS Solutions:

World's-first innovations available to those with single-sided deafness



## Augmented Xperience (AX) Platform

### Speech clarity without compromise

The Augmented Xperience platform introduces our new ear-to-ear (e2e) Wireless 4.0 which makes it easier for the brain to process speech sounds ensuring you hear clearer speech and a smoother sound experience. CROS compatibility now available in the level 2AX and 1AX.



Patients with single sided deafness can benefit from the AX technology with both rechargeable or traditional battery RIC options.



#### CROS Pure Charge&Go AX

Rechargeable RIC that delivers up to 29 hours of wear time



#### CROS Pure 312 AX

Traditional battery operated RIC that delivers up to 4 days of charge



#### CROS Styletto AX

Slim RIC that delivers up to 20 hours of wear time, including 5 hours of streaming



#### Features included in all technology levels:

- IP68 rating
- Directional pattern skewed toward the unaidable ear preserving CROS benefit even in noise
- Motion sensors provide improved speech understanding even while moving
- CROS mixing-level adjustment available in the Signia app for better balance
- Rocker switch for onboard controls AND app control of CROS settings\*\*
- Basic Remote Tuning through the TeleCare portal

#### Features included in premium technology levels:

- NEW Upgraded eWindScreen for smarter, faster wind noise reduction — levels 7/5/3/2 only\*\*\*\*
- NEW Improved AX Soundscape Processing to help you engage in challenging sound environments\*\*\* — levels 7/5/3 only
- Auto EchoShield to eliminate sound reverberation — level 7 only
- Own Voice Processing 2.0 (OVP2.0) in BiCROS fittings — levels 7/5/3 only
- HandsFree with CallControl for iOS devices — levels 7/5/3 only
- My WellBeing via the Signia app — levels 7/5/3 only

## Xperience Platform

### Ready to wear design

Signia delivers three CROS and BiCROS solutions on the Xperience platform.†

The revolutionary Signia Xperience platform offers the only ready-to-wear CIC, delivering clear speech understanding in an ultra-discreet size with instant-fit.



#### CROS Silk X

Traditional battery operated ready-to-wear CIC with up to 55 hours of wear time



#### Features include:

- Natural hearing experience and exceptional sound quality even while using headphones
- Ultra-discreet housing with click-sleeve options for individualized comfort
- Compatible with the Signia app and miniPocket™

#### LEGEND



\* Visit [signia-pro.com/en-us/claims-data](http://signia-pro.com/en-us/claims-data) for more details.

\*\* Onboard control not available on Styletto AX.

\*\*\* Upgraded eWindScreen requires Connexx 9.9 or later.

\*\*\*\* Available in products where OVP 2.0 has been activated. Requires Connexx 9.9 or later.

† Also available on the Xperience Platform: CROS Pure Charge&Go X, CROS Pure 312 X.

†† Signia Assistant is not available with CROS & BiCROS fittings