

# Powered by the Signia Augmented Xperience Platform



## Expanding on our mission to enhance human performance

Signia's latest innovations – the Pure 312 AX and the CROS Pure 312 AX – are powered by the new and cutting-edge Augmented Xperience (AX) platform, which delivers up to 96 channels with two dedicated processors for split processing and a host of clinically proven benefits, including new and enhanced audiological features. The platform is designed to deliver outstanding speech quality in quiet, moderate, and noisy environments!

## Provider Benefits

Fully featured RIC designed to **support a broad range of patients' needs:**

- Easily connect to iPhone® and ASHA compatible Android™ devices
- Integrated, compact RIC design with CROS option
- Rocker switch for onboard control
- Sleek design combined with new EarWear 3.0, to deliver higher physical fit rates and better retention for increased patient satisfaction
- Pure 312 AX's Zinc-Air battery delivers up to 4 days of wear time

## Wearer Benefits

- Sleek hearing aid design that is both comfortable and discreet (via EarWear 3.0)
- Fully featured hearing aid with optional support of induction looping systems and a CROS transmitter for single-sided deafness
- **Signia app compatibility** with access to Signia Assistant and more for on-demand, easy-to-use control
- **Pure 312 AX batteries** keep hearing aids running for up to four days before needing to be replaced
- **Own Voice Processing (OVP) 2.0** separates the wearer's own voice independently from other sounds - to increase own-voice comfort for easier communicating. **Plus! Improved AX Soundscape Processing** leverages the power of OVP 2.0 and uses own voice detection as another input to optimize how conversation partners are heard\*\*.
- **HandsFree calls**, via Bluetooth®, allows the wearer's iPhone to connect directly to their hearing aids for on-the-go conversations
- **Auto EchoShield** creates a cleaner, more immersive environment by neutralizing reverberation and sudden loud noises instantly, before they become an issue in atriums, hallways, lecture halls, and other complex environments
- **Upgraded eWindScreen†** reduces annoying wind sounds to allow wearers to better engage in outdoor communication or other outdoor activities by reacting smarter and faster to fluctuating wind noise.

## Key Features

- Outstanding speech clarity
- Android & iOS connectivity
- Signia Assistant available 24/7 through the Signia app



## Available Colors:



Signia app is equipped with the Signia Assistant for 24/7 support plus Signia Face Mask Mode Improves audibility for those communicating with someone wearing a face mask — all with the touch of a button in the Signia app.

\*Visit signia-pro.com/en-us/claims-data for more details.

\*\*Available in levels 3-7 when OVP 2.0 has been activated. Requires Connexx 9.9 or later.

\*\*\*Available in levels 2-7 in Connexx 9.9 or later.

**Type** **Pure 312 AX**

Matrix	S (110/46) M (119/60) P (122/65) HP (131/75)
--------	---

Battery Size	312 Zinc-Air
--------------	--------------

**Features**

Frequency channels/ compression controls	96/20 (48 channels per processor)
---	-----------------------------------

Hearing Programs	6
------------------	---

Bluetooth® Streaming	•
----------------------	---

TwinPhone <sup>1</sup>	•
------------------------	---

TeleCare™	•
-----------	---

Tinnitus therapies	•
--------------------	---

CROS/BiCROS <sup>2</sup>	•
--------------------------	---

Rocker Switch <sup>3</sup>	•
----------------------------	---

**App/Accessory Compatibility**

Signia App	•
------------	---

StreamLine TV/StreamLine Mic	•
------------------------------	---

miniPocket	•
------------	---

Noahlink Wireless	•
-------------------	---

Fitting set	EarWear 3.0
-------------	-------------

**Pure 312 AX works exclusively with the new EarWear 3.0.**

**Compatible with:**



StreamLine™ TV



StreamLine™ Mic



miniPocket™

<sup>1</sup> Activate via rocker switch, app or accessory

<sup>2</sup> Compatible with CROS Pure 312 AX

<sup>3</sup> Programmable via Connexx

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.



Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play and the Google Play logo are trademarks of Google Inc.

- Standard
- Optional

Signia is a registered trademark of WSAUD A/S. 04/23 SIV-15169-22 (16576)

