

Powered by the Signia Augmented Xperience Platform



Styletto AX



The Sound of Confidence

Signia's latest innovation – the Styletto AX with CROS compatibility – is powered by the cutting-edge Augmented Xperience (AX) platform. This award-winning design is preferred by 8 out of 10 hearing aid wearers.*

Key Sound Benefits

- **Auto EchoShield** creates a cleaner, more immersive environment by neutralizing reverberation and sudden loud noises instantly, before they become an issue.
- **Own Voice Processing (OVP™) 2.0** utilizes split processing of the Augmented Xperience platform to separate the wearer's own voice from simultaneous ambient sounds. **Plus! Improved AX Soundscape Processing** leverages the power of OVP 2.0 and uses own voice detection as another input to optimize how conversation partners are heard.**
- **HandsFree for iOS** allows wearers to use their hearing aids for hands-free calls with their iPhone® in different environments.
- Powered by the **first-of-its-kind AX platform** with two dedicated processors.
 - **More than 25% better speech understanding** in noise compared to own hearing aids¹.
 - **94% of wearers report excellent speech understanding** with Signia AX² New Features.
- **Upgraded eWindScreen***** reduces annoying wind sounds to allow wearers to better engage in outdoor communication or other outdoor activities by reacting smarter and faster to fluctuating wind noise.

Provider Benefits

- Its refined **slim line hearwear™** styling enables you to attract new and existing hearing aid wearers alike.
- Featured **rechargeable SLIM-RIC** models are designed to support a broad range of patients' needs.
- **Innovative SLIM-RIC design** for an even more comfortable fit.
- CROS option available.
- When presented with a choice, **84%** of the survey participants reported they would visit the office that included Styletto in their display versus those who did not.*
- Attract new patients with a stigma-fighting design.
- **TeleCare™** ready.

Wearer Benefits

- **Innovative SLIM-RIC design** for an even more comfortable, discreet fit.
- Easily connects to iPhone and ASHA-compatible Android™ devices for effortless streaming from any Bluetooth device.
- **Wireless charging** with Qi technology offers ease of use.
- Up to **20 hours** of runtime on a single charge, or 17 hours with 5 hours of streaming, plus enjoy a 30-minute fast charge for an additional 5.5 hours of hearing aid power.
- **Portable**, wireless charging case that holds up to 4 additional days of charges.

Available Colors:



Signia app is equipped with the Signia Assistant for 24/7 support plus Signia My WellBeing. My WellBeing provides health insights to track steps, movement, wear time and more to support your patients overall health.

* Hakvoort, Burton: Increasing Style, Reducing Stigma: The Styletto Solution (Signia Backgrounder, 2018)

** Available in levels 3-7 when OVP 2.0 has been activated. Requires Connexx 9.9 or later.

*** Available in levels 2-7 in Connexx 9.9 or later.

¹ Jensen N.S., Pischel C., Taylor B. & Schulte M. 2021b. Performance of Signia AX in at-home listening situations. Signia White Paper.

² Visit signia-pro.com/en-us/claims-data for more details

Type	Styletto AX
Available Performance Level	7
Matrix	S (110/46) M (119/60) P (122/65)
Battery Size	Lithium-ion power cell
Features	
Frequency channels/ compression controls	96/20 (48 channels per processor)
Hearing Programs	6
Bluetooth® Streaming	•
TwinPhone ¹	•
TeleCare	•
Tinnitus therapies	•
CROS/BiCROS ²	•
Rechargeable	•
HandsFree with CallControl for iOS	•
App/Accessory Compatibility	
Signia app	•
StreamLine™ TV/StreamLine™ Mic	•
miniPocket™	•
Noahlink Wireless	•
Fitting set	EarWear 3.0

Styletto AX & CROS Styletto AX Charging for Easy Operation



Default

Styletto AX
Charging Case



Styletto AX
Dry&Clean Charger

Compatible with:



StreamLine™ TV



StreamLine™ Mic



miniPocket™

¹ Activate via app or accessory

² Compatible with CROS Styletto AX

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

- Standard



Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play and the Google Play logo are trademarks of Google Inc.

**Be
Brilliant™**