

Powered by the Signia Augmented Xperience Platform



Insignia Charge&Go AX



Expanding on our mission of enhancing human performance

The Insignia Charge&Go AX is powered by the world's first Augmented Xperience platform which delivers up to 96 channels with two dedicated processors for split processing and a host of clinically proven benefits including NEW and enhanced audiological features. The Augmented Xperience platform is designed to deliver outstanding speech quality in quiet, moderate, and noisy environments!

Provider Benefits

- Featured rechargeable ITC and ITE models designed to support a broad range of patient's needs:
 - **Easily connect to iPhone® and ASHA-compatible Android™ devices**
 - **Convenient push-button for onboard control**
- Small design for a higher fitting rate and all-day comfort for increased patient satisfaction
- Easy insertion and removal to address challenges with dexterity
- First **contactless charger for easy operation**
 - No specific charging position required
 - Wearers intuitively place the Insignia Charge&Go hearing aids in the charger

Wearer Benefits

- **Small, custom-made hearing aids tailored to the wearer's ear canal providing all-day wearing comfort**
- **Easily connects to iPhone and ASHA compatible Android devices** for effortless streaming from any Bluetooth® device
- Contactless charger offers ease of use
- Innovative charger provides up to 24 hours of use on a single charge
- **Hands-free calls**, via Bluetooth®, allows the wearer's iPhone to connect directly to their hearing aids for clearer conversations
- **Auto EchoShield** creates a cleaner, more immersive environment by neutralizing reverberation and sudden loud noises instantly, before they become an issue in atriums, hallways, lecture halls, and other complex environments
- **Upgraded eWindScreen***** reduces annoying wind sounds to allow wearers to better engage in outdoor communication or other outdoor activities by reacting smarter and faster to fluctuating wind noise.

Key Features*

- Built in Li-ion Cell for contactless charging
- Powered by the **World's First Augmented Xperience platform with two dedicated processors**
 - More than **25% better speech understanding in noise** compared to own hearing aids**
 - **94%** of wearers report **excellent speech understanding** with **Signia AX****
- TeleCare™ ready



Custom Design



iOS & Android Connectivity



Contactless Charging



Tinnitus Relief



Automatic Situation Detection



IP 68



Signia App



24/7 Signia Assistant



HandsFree with CallControl for iOS



Auto EchoShield



My WellBeing



Upgraded eWindScreen

Available Colors:



Black



Dark Brown



Brown



Beige American Flag



Tan



Beige



Mocha



Black American Flag

Shells available in all six colors with clear, red, and blue available by request.

Signia app is equipped with Signia Face Mask Mode.

It improves audibility for those communicating with someone wearing a face mask — all with the touch of a button in the Signia app.

* Refer to Feature Overview Chart in the AX Platform Brochure

** Visit signia-pro.com/en-us/claims-data for more details

Type	Insio Charge&Go AX
Available Performance Level	7
Matrix	124/65 118/55 113/50
Battery Size	Lithium-ion power cell
Features	
Frequency Controls/ Compression Channels	7AX - 20/96 (48 per processor)
Hearing Programs	6
Bluetooth® Streaming	•
TwinPhone ¹	•
TeleCare	•
Tinnitus therapies	•
Push Button	•
Rechargeable	•
App/Accessory Compatibility	
Signia App	•
StreamLine TV/StreamLine Mic	•
miniPocket	•
Noahlink Wireless	•

Insio Charge&Go AX with the world's first contactless charger for easy operation



Insio Charger

- Innovative charger fully charges in only 4 hours and provides up to 24 hours of use on a single charge
- No specific charging position required
- Wearers intuitively place the Insio Charge&Go hearing aids in the charger

Compatible with:



StreamLine™ TV



StreamLine™ Mic



miniPocket™

¹ Activate via push button, app or accessory

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.



Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play and the Google Play logo are trademarks of Google Inc.