

Powered by the Signia Augmented Xperience Platform



Delivering on our mission of enhancing human performance

The Pure Charge&Go T AX and the CROS Pure Charge&Go T AX are powered by the new and cutting-edge Augmented Xperience (AX) platform, which delivers up to 96 channels with two dedicated processors for split processing and a host of clinically proven benefits, including new and enhanced audiological features. The platform is designed to deliver outstanding speech quality in quiet, moderate, and noisy environments!

Provider Benefits

- Fully featured RIC designed to **support a broad range of patient's needs:**
 - Easily connect to iPhone® and ASHA compatible Android™ devices**
 - Integrated **telecoil in a compact RIC design and CROS options**
 - Rocker switch for onboard control
 - Sleek design combined with new EarWear 3.0 which is proven to deliver higher physical fit rates and better retention for increased patient satisfaction
- Rechargeable with **multiple charger options** to choose from to appeal to each patient's unique needs
 - Portable Charger** – provides **3 days of extra charges** for patients on the go
 - Pure Charger** – simplifies charging with its **small, freestanding design**
 - Dry&Clean Charger** – to **keep hearing aids in top working conditions**

Wearer Benefits

- Sleek hearing aid design** that is both comfortable and discreet (via EarWear 3.0)
- Multiple charger options to match your lifestyle with **up to 36 hours of wear time on a single charge.**
- HandsFree calls**, via Bluetooth®, allows the wearer's iPhone to connect directly to their hearing aids for clearer conversations
- Own Voice Processing (OVP) 2.0** separates the wearer's own voice independently from other sounds - to increase own-voice comfort for easier communicating. **Plus! Improved AX Soundscape Processing** leverages the power of OVP 2.0 and uses own voice detection as another input to optimize how conversation partners are heard.†
- Auto EchoShield** creates a cleaner, more immersive environment by neutralizing reverberation and sudden loud noises instantly, before they become an issue in atriums, hallways, lecture halls, and other complex environments
- Upgraded eWindScreen**†† reduces annoying wind sounds to allow wearers to better engage in outdoor communication or other outdoor activities by reacting smarter and faster to fluctuating wind noise.

Key Features*

- Up to 36 hours of use on a single charge with telecoil, and up to 24 hours without telecoil (reflects 5 hours of streaming in both options)
- Powered by the **world's first Augmented Xperience platform with two dedicated processors**
 - More than **25% better speech understanding in noise** compared to own hearing aids**
 - 94%** of wearers report **excellent speech understanding with Signia AX****
- TeleCare™ Ready

HandsFree with CallControl for iOS	Bluetooth Connectivity	Li-ion Rechargeability	Upgraded eWindScreen††
Own Voice Processing 2.0	Dynamic Soundscape Processing 2.0	Acoustic Motion Sensors	Improved AX Soundscape Processing†
Auto EchoShield	Tinnitus Relief	Telecoil	Signia App
	24/7 Signia Assistant	My Well Being	

Available Colors:



Signia app is equipped with the Signia Assistant for 24/7 support plus Signia Face Mask Mode Improves audibility for those communicating with someone wearing a face mask — all with the touch of a button in the Signia app.

* Refer to Feature Overview Chart in the AX Platform Brochure

** Visit signia-pro.com/en-us/claims-data for more details.

† Available in levels 3-7 when OVP 2.0 has been activated. Requires Connexx 9.9 or later.

†† Available in levels 2-7 in Connexx 9.9 or later.

Type	Pure Charge&Go T AX
Available Performance Level	7
Matrix	S (110/46) M (119/60) P (122/65) HP (131/75)
Battery Size	Lithium-ion power cell
Features	
Frequency channels/ compression controls	96/20 (48 channels per processor)
Hearing Programs	6
Bluetooth® Streaming	•
TwinPhone ¹	•
TeleCare™	•
Tinnitus therapies	•
CROS/BiCROS ²	•
Rocker Switch ³	•
Rechargeable	•
Telecoil	◦
App/Accessory Compatibility	
Signia App	•
StreamLine TV/StreamLine Mic	•
miniPocket	•
Noahlink Wireless	•
Fitting set	EarWear 3.0

Signia is a registered trademark of WSAUD A/S. SIV-15179-22 (16604)

Pure Charge&Go T AX works exclusively with the new EarWear 3.0.

Multiple charger options available



Pure Portable Charger
Standard Option
Default



Pure Charger
Standard Option



Pure Dry&Clean Charger
Standard Option

Compatible with:



StreamLine™ TV



StreamLine™ Mic



miniPocket™

¹ Activate via rocker switch, app or accessory

² Compatible with CROS Pure Charge&Go T AX

³ Programmable via Connexx

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

- Standard
- Optional



Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Android, Google Play and the Google Play logo are trademarks of Google Inc.

**Be
Brilliant™**