



Styletto AX Charger Styletto IX Charger

Data Sheet



Styletto AX / IX Charger | Technical Data

Type

AX

IX

Device property

Dimensions (LxWxH)	70.8 x 49.9 x 29.5 mm
Weight without cable	50 g
Ingress Protection Rating	IP XO (ordinary)
Protection against electric shock	class III

Electrical property – Charger

Power supply connector	USB-C
Input voltage	5 V DC
Input current	≤ 0.5 A

Electrical property – Power supply adapter

Input voltage	100 – 240 V AC
Input current	0.15 – 0.3 A
Input frequency	50 – 60 Hz
Output voltage	5 V DC
Output current	1 A

Built-in battery

Battery type	Lithium-ion	
Battery capacity (see label on device)	400 mAh / 450 mAh	450 mAh
Battery runtime	4 days in typical use	
Charging time (USB)	~ 3 h	
Charging time (wireless)	~ 6 h	

Operating conditions

Temperature (USB charging)	+10 ... +35 °C / +50 ... +95 °F
Temperature (wireless charging)	+10 ... +30 °C / +50 ... +86 °F
Relative humidity	20 ... 75 %
Minimum air pressure	700 hPa

Storage conditions

Temperature (recommended)	+10 ... +25 °C / +50 ... +77 °F
Temperature (max.)	+10 ... +40 °C / +50 ... +104 °F
Relative humidity (recommended)	20 ... 80 %
Relative humidity (max.)	10 ... 80 %

Transportation conditions

Temperature	- 20 ... +60 °C / - 4 ... +140 °F
Relative humidity	5 ... 90 %

Note: Due to the design of the Styletto AX / IX, not all earwear couplings will fit inside the Styletto AX / IX Charger.

Recommended couplings: S, M and P-Receiver (length 0 – 3) with Earware 3.0 open and tulip (all sizes) or sleeves (size XS – M).

Styletto AX / IX Charger | Conformance information

		Safety	EN 62368-1:2020
EU	CE	EMC	EN 60601-1-2:2015 + A1:2021
USA	FCC		47 CFR Part 15
Canada	IC		ICES-003 Issue 7 / 2020
Korea	KCC		KN 32 / KN 35