

GOVERNMENT SERVICES RIC IX/AX ORDER FORM (MAY 2024)

TEL: (800) 955-4327 FAX: (732) 562-6653 EMAIL: va.customercare@signiausa.com



Account #

Date:

Bill To Address:

Contact:
Tel:

Ship To Address:

Email:

P.O. No.

Patient Information Last Name, First Name, Middle Initial Last 4 S.S.#

Shipping: No charge UPS 2nd Day
Warranty: 36 months

Hearing Thresholds		250 Hz	500 Hz	750 Hz	1 kHz	1.5 kHz	2 kHz	3 kHz	4 kHz	6 kHz	8 kHz
Right											
Left											

Earmolds will ship with the RIC when ordered together

Check if RIC Ordered with Earmolds

1 Model	Pure Charge&Go TIX	CROS Pure Charge&Go IX	Styletto IX	CROS Styletto IX	CROS Pure Charge&Go AX	Pure 312 AX	CROS Pure 312 AX	Styletto AX	CROS Styletto AX
	R L	R L	R L	R L	R L	R L	R L	R L	R L
Hearing instrument	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>
Battery	Li-ion	Li-ion	Li-ion		Li-ion	ZA 312	ZA 312	Li-ion	
Program selector	Rocker	Rocker	via Signia App or miniPocket Remote		Rocker	Rocker	Rocker	via Signia App or miniPocket Remote	
Volume control	Rocker	Rocker	via Signia App or miniPocket Remote		Rocker	Rocker	Rocker	via Signia App or miniPocket Remote	

2 Receiver 3.0 Power Level

	R	L
S - 110/46 dB	<input type="checkbox"/>	<input type="checkbox"/>
M - 119/60 dB	<input type="checkbox"/>	<input type="checkbox"/>
P - 122/65 dB	<input type="checkbox"/>	<input type="checkbox"/>
HP* - 131/75 dB	<input type="checkbox"/>	<input type="checkbox"/>

*Not available for Styletto IX/AX, CROS Styletto IX/AX
*Custom impression required

3 Receiver Lengths

	R	L
0	<input type="checkbox"/>	<input type="checkbox"/>
1	<input type="checkbox"/>	<input type="checkbox"/>
2	<input type="checkbox"/>	<input type="checkbox"/>
3	<input type="checkbox"/>	<input type="checkbox"/>
4*	<input type="checkbox"/>	<input type="checkbox"/>
5*	<input type="checkbox"/>	<input type="checkbox"/>

*Not recommended for Styletto IX/AX

6 Accessories
(Additional charge for accessories)

StreamLine™ TV	<input type="checkbox"/>
StreamLine™ Mic	<input type="checkbox"/>
miniPocket™	<input type="checkbox"/>



4 Colors

Pure RICs	R	L	Styletto (first color denotes hearing aid color; second color denotes spine color)	R	L
Beige (1)	<input type="checkbox"/>	<input type="checkbox"/>	Black/Black	<input type="checkbox"/>	<input type="checkbox"/>
Black (5)	<input type="checkbox"/>	<input type="checkbox"/>	Beige/Beige	<input type="checkbox"/>	<input type="checkbox"/>
Dark Champagne (37)	<input type="checkbox"/>	<input type="checkbox"/>	Deep Brown/ Deep Brown	<input type="checkbox"/>	<input type="checkbox"/>
Deep Brown (72)	<input type="checkbox"/>	<input type="checkbox"/>	Black/Fine Gold*	<input type="checkbox"/>	<input type="checkbox"/>
Fine Gold (647)	<input type="checkbox"/>	<input type="checkbox"/>	Black/Graphite*	<input type="checkbox"/>	<input type="checkbox"/>
Graphite (646)	<input type="checkbox"/>	<input type="checkbox"/>	Black/Silver	<input type="checkbox"/>	<input type="checkbox"/>
Pearl White (7)	<input type="checkbox"/>	<input type="checkbox"/>	Cosmic Blue/Fine Gold*	<input type="checkbox"/>	<input type="checkbox"/>
Rose Gold (927)	<input type="checkbox"/>	<input type="checkbox"/>	Cosmic Blue/Rose Gold	<input type="checkbox"/>	<input type="checkbox"/>
Sandy Brown (31)	<input type="checkbox"/>	<input type="checkbox"/>	Cosmic Blue/Silver	<input type="checkbox"/>	<input type="checkbox"/>
Silver (926)	<input type="checkbox"/>	<input type="checkbox"/>	Snow White/Fine Gold*	<input type="checkbox"/>	<input type="checkbox"/>
			Snow White/Rose Gold	<input type="checkbox"/>	<input type="checkbox"/>
			Snow White/Silver	<input type="checkbox"/>	<input type="checkbox"/>
			Forest Camo/Black	<input type="checkbox"/>	<input type="checkbox"/>
			Desert Camo/Black	<input type="checkbox"/>	<input type="checkbox"/>
			Snow White/Snow White Gloss**	<input type="checkbox"/>	<input type="checkbox"/>

*N/A on Styletto IX **For Styletto IX only

7 Chargers

Pure (Portable Charger is default)	Styletto (Styletto Charger is default)
Portable Pure Charger <input type="checkbox"/>	Styletto Charger <input type="checkbox"/>
Dry&Clean Charger <input type="checkbox"/>	Styletto Dry and Clean Charger <input type="checkbox"/>
Pure Standard Charger <input type="checkbox"/>	

8 Fitting Formula

NAL-NL2 Fitting Formula	<input type="checkbox"/>
IX Fitting Formula	<input type="checkbox"/>
AX Fitting Formula	<input type="checkbox"/>

5 Ear Tips and Sleeves

	R	L
Ear Tips 3.0		
Open S(5 mm)	<input type="checkbox"/>	<input type="checkbox"/>
Open M(7 mm)	<input type="checkbox"/>	<input type="checkbox"/>
Open L(10 mm)	<input type="checkbox"/>	<input type="checkbox"/>
Tulip S	<input type="checkbox"/>	<input type="checkbox"/>
Tulip M	<input type="checkbox"/>	<input type="checkbox"/>
Round Dome 1 S	<input type="checkbox"/>	<input type="checkbox"/>
Round Dome 1 M	<input type="checkbox"/>	<input type="checkbox"/>
Round Dome 1 L	<input type="checkbox"/>	<input type="checkbox"/>
Sleeves 3.0		
Vented XS	<input type="checkbox"/>	<input type="checkbox"/>
Vented S	<input type="checkbox"/>	<input type="checkbox"/>
Vented M	<input type="checkbox"/>	<input type="checkbox"/>
Vented L	<input type="checkbox"/>	<input type="checkbox"/>
Closed XS	<input type="checkbox"/>	<input type="checkbox"/>
Closed S	<input type="checkbox"/>	<input type="checkbox"/>
Closed M	<input type="checkbox"/>	<input type="checkbox"/>
Closed L	<input type="checkbox"/>	<input type="checkbox"/>
Power XS	<input type="checkbox"/>	<input type="checkbox"/>
Power S	<input type="checkbox"/>	<input type="checkbox"/>
Power M	<input type="checkbox"/>	<input type="checkbox"/>
Power L	<input type="checkbox"/>	<input type="checkbox"/>

Special Instructions:

Please Include a Paper User Manual with these hearing aids

GOVERNMENT SERVICES RIC ORDER FORM (MAY 2024)

TEL: (800) 955-4327 FAX: (732) 562-6653 EMAIL: va.customercare@signiausa.com



RIC Earmold 3.0 and RIC Embedded Earmold 3.0			
RIC Earmold 3.0			
	R L		
Earmold Style RIC Earmold 3.0			
Hollow			
RIC	<input type="checkbox"/>	<input type="checkbox"/>	
RIC Lock	<input type="checkbox"/>	<input type="checkbox"/>	
Half Shell	<input type="checkbox"/>	<input type="checkbox"/>	
Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	
Solid			
RIC	<input type="checkbox"/>	<input type="checkbox"/>	
RIC Lock	<input type="checkbox"/>	<input type="checkbox"/>	
Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	
Full Shell	<input type="checkbox"/>	<input type="checkbox"/>	
Half Shell	<input type="checkbox"/>	<input type="checkbox"/>	
Semi-Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	
Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	
RIC Embedded Earmold 3.0 (Only available in RIC, RIC Lock, and Skeleton Lock)			
	R L		
Earmold Style RIC Embedded Earmold 3.0			
RIC	<input type="checkbox"/>	<input type="checkbox"/>	
RIC Lock	<input type="checkbox"/>	<input type="checkbox"/>	
Skeleton Lock	<input type="checkbox"/>	<input type="checkbox"/>	
Side Indicators (pick only 1 type)		RIC Earmold 3.0	RIC Embedded Earmold 3.0
Large R/L Letters	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Large Dot	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Patient Name	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Removal String		RIC Earmold 3.0	RIC Embedded Earmold 3.0
Removal String	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Wax Options		RIC Earmold 3.0	RIC Embedded Earmold 3.0
Nanocare	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Defaults are in bold.			

Material	RIC Earmold 3.0	RIC Embedded Earmold 3.0
	R L	R L
Hard	<input type="checkbox"/>	<input type="checkbox"/>
Soft (N/A for Hollow and Embedded Earmolds)	<input type="checkbox"/>	<input type="checkbox"/>
Matte or Gloss finish	RIC Earmold 3.0	RIC Embedded Earmold 3.0
Matte <input type="checkbox"/>		Gloss <input type="checkbox"/>
Clear	<input type="checkbox"/>	<input type="checkbox"/>
Rose	<input type="checkbox"/>	<input type="checkbox"/>
Beige	<input type="checkbox"/>	<input type="checkbox"/>
Brown (hard material only)	<input type="checkbox"/>	<input type="checkbox"/>
Medium Brown	<input type="checkbox"/>	<input type="checkbox"/>
Dark Brown	<input type="checkbox"/>	<input type="checkbox"/>
Silicone with gloss finish	RIC Earmold 3.0	RIC Embedded Earmold 3.0
Green	<input type="checkbox"/>	<input type="checkbox"/>
Pink	<input type="checkbox"/>	<input type="checkbox"/>
Clear Glitter	<input type="checkbox"/>	<input type="checkbox"/>
Purple Glitter	<input type="checkbox"/>	<input type="checkbox"/>
Red	<input type="checkbox"/>	<input type="checkbox"/>
Orange	<input type="checkbox"/>	<input type="checkbox"/>
Yellow	<input type="checkbox"/>	<input type="checkbox"/>
Blue	<input type="checkbox"/>	<input type="checkbox"/>
Lilac	<input type="checkbox"/>	<input type="checkbox"/>
Black	<input type="checkbox"/>	<input type="checkbox"/>
White	<input type="checkbox"/>	<input type="checkbox"/>
Canal Length	RIC Earmold 3.0	RIC Embedded Earmold 3.0
Short	<input type="checkbox"/>	<input type="checkbox"/>
Medium	<input type="checkbox"/>	<input type="checkbox"/>
Long (hard solid material only)	<input type="checkbox"/>	<input type="checkbox"/>
Venting	RIC Earmold 3.0	RIC Embedded Earmold 3.0
Open (Hollow only)	<input type="checkbox"/>	<input type="checkbox"/>
No Vent	<input type="checkbox"/>	<input type="checkbox"/>
Pressure Vent 1.0mm	<input type="checkbox"/>	<input type="checkbox"/>
1.2 mm	<input type="checkbox"/>	<input type="checkbox"/>
1.6 mm	<input type="checkbox"/>	<input type="checkbox"/>
SAV 2.0 mm	<input type="checkbox"/>	<input type="checkbox"/>
SAV 3.0 mm	<input type="checkbox"/>	<input type="checkbox"/>
As large as possible	<input type="checkbox"/>	<input type="checkbox"/>
Trench Vent (hard material only)	<input type="checkbox"/>	<input type="checkbox"/>
Semi IROS	<input type="checkbox"/>	<input type="checkbox"/>

HP 3.0 Earmold (Not available with Styletto IX or AX)	
HP Earmold (131/75 dB)	
Material (Only available in hard material)	
Faceplate Colors	
Beige	R L
Mocha	<input type="checkbox"/>
Brown	<input type="checkbox"/>
Shell Colors	
Clear	<input type="checkbox"/>
Tan	<input type="checkbox"/>
Beige	<input type="checkbox"/>
Brown	<input type="checkbox"/>
Red	<input type="checkbox"/>
Blue	<input type="checkbox"/>
Venting (HP Only)	
No vent	<input type="checkbox"/>
Pressure Vent 1.0mm	<input type="checkbox"/>
1.6mm	<input type="checkbox"/>
Shell Options	
Retrieval Line	<input type="checkbox"/>
Canal Lock	<input type="checkbox"/>
Helix Lock	<input type="checkbox"/>
Skeleton Lock	<input type="checkbox"/>
Canal Length	
Medium	<input type="checkbox"/>
Long	<input type="checkbox"/>
Deep	<input type="checkbox"/>
Side Indicators (pick only 1 type)	
Large R/L Letters	<input type="checkbox"/>
Large Dot	<input type="checkbox"/>
Patient Name	<input type="checkbox"/>
Wax Protection (HP Only)	
Nanocare Cerustop XL	<input type="checkbox"/>

Copyright © 2024 Signia is a registered trademark of WSAUD A/S. Part #10536367 16.0 3/24 SIV-18648-24



For internal use only

Aids: 1 2 Earmolds: 1 2
 Receiver: 1 2 Impressions: 1 2
 Charger: 1 2 Remotes: 1 2

Impression: Yellow Blue Pink Peach White Purple Trans Brown Turq Green Org

RECEIVED = NJ Bag # _____ ID # _____