

About Signia

Only Signia combines industry-leading audiology and technology to deliver iconic innovations in hearing to enhance human performance. We put hearing aid wearers first, with products that fit their wants, needs and style. The wearer is at the center of our universe and their needs are woven into our DNA.

For us, this is not about correcting a loss, it's about gaining an edge. Not just performing, but performing at their best.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by WSAUD A/S is under license. Other trademarks and trade names are those of their respective owners. Other trademarks and trade names are those of their respective owners.

Legal Manufacturer

Sivantos, Inc.
10 Constitution Avenue
Piscataway, NJ 08854
Phone: (800) 766-4500
Fax: (732) 562-6696



Copyright © 2022 Signia is a registered trademark of WSAUD A/S.
4/22 SIV-12982-22

Be
Brilliant™

* Visit signia-pro.com/en-us/claims-data for more details.

Find us on:



signia.net/en-us



facebook.com/signiahearingUS



youtube.com/signiahearing



instagram.com/signiahearingus

Be
Brilliant™



Immerse yourself
in outstanding
speech clarity

Signia Augmented Xperience (AX)



signia

Pure Charge&Go T AX
Pure 312 AX
Insignia Charge&Go AX





Speech clarity without compromise

Our range of Signia Augmented Xperience (AX) hearing aids change the way you hear the world.

No more straining to hear speech above background noise: Signia's revolutionary AX technology lets you follow conversations with ease, even in difficult listening situations. It helps you stay focused on the conversation you want to hear, while placing background noises where they belong – in the background.

Enjoy Hearing with Greater Ease



3 out of 4 wearers experience less listening effort with AX than with their own hearing aids*



Augmented Xperience



Immersive Soundscape

Revolutionizing the way you process sound

Traditionally, hearing aids have processed all sounds under one set of rules, resulting in signals competing against each other. Now, Signia Augmented Xperience (AX) hearing aids can separate speech from background noise using Augmented Focus™ technology to bring the sounds you want to hear into focus.

Two dedicated processors can fully optimize different sounds without compromising one or the other. The sounds are then recombined to deliver outstanding speech clarity in a fully immersive soundscape.


This enhanced technology helps you with more than just hearing. 7 out of 10 people say they feel more confident with AX than with their previous hearing aids, and 80% feel more energized!*

Excellent speech understanding with Signia AX



94% of wearers report excellent speech understanding with Signia AX*




Automatic
Situation Detection


e2e
wireless 4.0


Tinnitus Relief

AX benefits at a glance

Automatic Situation Detection

Advanced sensor technology recognizes your movement and automatically adjusts your settings for better understanding of speech from any direction and maximizes your audibility as you move throughout the day.

Tinnitus Relief


Signia AX features three different tinnitus relief strategies, including static noise therapy, ocean wave therapy, and exclusive Notch therapy, to find the right solution to reduce or remove annoying ringing or buzzing in your ears.

e2e wireless 4.0

A great leap for true binaural hearing through a synchronization of sound between ears to ensure clearer speech processing with improved awareness of where sounds are coming from. Plus, e2e wireless 4.0 enhances feedback cancellation to avoid hearing aid squealing.



Better Than Normal Hearing*




OVP 2.0


HandsFree for iOS


Auto EchoShield

AX is enhancing your hearing experience with revolutionary new advancements

Own Voice Processing (OVP) 2.0**

The newly developed Own Voice Processing 2.0 utilizes split processing to help you hear yourself talking more naturally.

HandsFree for iOS

Use your hearing aids for hands-free calls with select smartphones in different environments. The hearing aids' built-in microphones let your conversation partners hear your voice in high quality.

Auto EchoShield

Auto EchoShield delivers automatic, real time echo reduction that ensures listening comfort in all situations, even in halls and rooms with large amounts of background noise.***

** Not available in Insio Charge&Go AX

*** Level 7 AX technology

Signia AX With Own Voice Processing 2.0
Ensures Great Sound of Own Voice



88%
of wearers like their own voice with
Signia AX with Own Voice Processing 2.0*



CROS



Rechargeable



Slim Top Design

Pure Charge&Go T AX

Signia's smallest rechargeable hearing aid with telecoil and a slimmer, sleek design.

- Portable rechargeability
- Bluetooth® streaming to iPhone® and Android™
- CROS option available
- Rocker switch for onboard control

Pure Charge&Go T AX is available in ten colors:



Black




Graphite

Dark
Champagne

Silver

Pearl
WhiteFine
GoldDeep
BrownSandy
BrownRose
Gold

Beige



CROS

Exchangeable
Batteries

Pure 312 AX

Signia's latest innovation, ideal for both experienced and first-time hearing aid wearers. Pure 312 AX offers a fresh, elegantly slim design, powered by long-lasting exchangeable batteries.

- Powered by replaceable battery, delivering a wear-time of up to 4 days
- Bluetooth® streaming to iPhone® and Android™
- CROS option available
- Rocker switch for onboard control

Pure 312 AX is available in ten colors:



Black




Graphite

Dark
Champagne

Silver

Pearl
WhiteFine
GoldDeep
BrownSandy
BrownRose
Gold

Beige



Customized Fit



Rechargeable

Inσιο Charge&Go AX

The world's first rechargeable custom hearing aids with contactless charging and Bluetooth® connectivity.

- Up to 24 hours of use on a single charge
- Bluetooth® streaming to iPhone® and Android™

Inσιο Charge&Go AX is available in six colors designed to match skin tone:



Black



Dark Brown



Brown



Mocha



Tan



Beige



iPhone & Android Connectivity



Recharge on-the-go

Connectivity meets convenience

All AX models are packed with streaming functionality to listen to music, take phone calls, and stream audio – from iPhone® and Android™ alike. Plus, enjoy handsfree calls with select iOS phones.

Convenient charging

AX is available in a variety of battery and rechargeable options, including Inσιο Charge&Go AX, the world's first custom hearing aids with contactless charging for ease of use.



Pure Charge&Go T AX Portable Charger



Pure Charge&Go T AX Dry&Clean Charger



Pure Charge&Go T AX Charger



Inσιο Charge&Go AX Charger

24/7 Signia
AssistantFace Mask
Mode

TeleCare

Personal support at your fingertips with Signia Assistant


Signia is committed to helping you live your life without interruptions – that means on-the-go support and leading technology all in one. Through Signia's advanced app, you have access to hearing aid controls, streaming capabilities with multiple devices, and so much more. Plus, you can easily connect to both iPhone® and Android™.

- On demand support via Artificial Intelligence Signia Assistant, available 24/7
- Improved speech intelligibility when talking with someone wearing a mask at the touch of a button with Face Mask Mode
- Virtual connection to your hearing healthcare provider for remote care services via Signia TeleCare

Based on the data of more than 50K users, we further improved the “brain” of the Signia Assistant machine learning to give users even better support.



93%
of hearing aid wearers report that Signia Assistant raises their satisfaction with hearing aids in difficult listening situations.*



Download the Signia App today to see the convenience and ease of use.



Compatible with iOS and Android. The app is free of charge.