
WAX GUARD

Quick Guide

Wax Guard Systems		
WG 3.0 NanoCare 10993649 	Use with: - miniReceiver 3.0 - ThinTube 3.0 and Power ThinTube 3.0 - Ear Mould 3.0	 *Note: All tubing/miniReceivers are fitted with a wax guard. For custom moulds an additional wax guard will be added if the canal length is long (as pictured).
WG 3.0 CeruSTOP XL 10999334 	Use with: - miniReceiver 3.0 HP Shell	

WAX GUARD

Quick Guide

Wax Guard Systems		
miniR Wax Guard 10962332 	Fits in the end of : - miniReceiver 2.0 - Silk - Active	
QuickGuard Pack 10962504 	Fits in the end of : - Click Mould 2.0 - HP Shell 2.0 - Insio Hearing Aids* *All models including Insio Charge&Go AX.	